


DIALYSIS SET


single lumen "Y"



Set components:

- single lumen catheter "Y"
- needle 18 G, length 7 cm
- guidewire
- dilator I
- dilator II

Optional components



Reference: guidewire "J" type	Reference: guidewire Straight type	Catheter size [F]	Catheter length [cm]	Dilator I/II [F]/[F]	Guidewire size [inch]	Guidewire length [cm]
ZDJY8F10	ZDJY8F10P	8	10	8/10	0.035"	60
ZDJY8F15	ZDJY8F15P		15			
ZDJY8F18	ZDJY8F18P		18			
ZDJY8F20	ZDJY8F20P		20			

When ordering use Ref. Number, e.g.: ZDJY8F20.

NOTE:

On demand, it is possible to select additional components.

While ordering, add the following endings to the reference index:

S – scalpel

Y – Y hub

SR – Guiding syringe 5 ml

SL – Syringe 5 ml

KLLT – TEGO connector

KLLI – Injection cap

IY – V needle

N – Nitinol guidewire

H – Hydrophilic coating

E.g.: ZDJY8F20HPYKLLT.

Different dialysis catheter sizes and lengths available on demand. For detailed information, please contact the company's representative.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.

