

DIALYSIS SET

triple lumen



Set components:

- triple lumen catheter
- needle 18 G, length 7 cm
- guidewire
- dilator I
- dilator II

Optional components



Straight type

| Reference: | Catheter size [F] | Catheter length [cm] | Dilator I/II [F]/[F] | Guidewire size [inch] | Guidewire length [cm] |
|------------|-------------------|----------------------|----------------------|-----------------------|-----------------------|
| ZDT12F15 | 12 | 15 | 10/12 | 0.035" | 70 |
| ZDT12F18 | | 18 | | | |
| ZDT12F20 | | 20 | | | |

When ordering use Ref. Number, e.g.: ZDT12F15.

NOTE:

On demand, it is possible to select additional components.

While ordering, add the following endings to the reference index:

S – scalpel

Y – Y hub

SR – Guiding syringe 5 ml

SL – Syringe 5 ml

KLLT – TEGO connector

KLLI – Injection cap

IY – V needle

N – Nitinol guidewire

H – Hydrophilic coating

E.g.: ZDT12F15HKLLIYN.

Different dialysis catheter sizes and lengths available on demand. For detailed information, please contact the company's representative.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.

