

## ONE-STEP METHOD DRAINAGE SET



- Set components:
- straight or pigtail type catheter
  - two pieces needle
  - band clip
  - collar

Straight or pigtail catheter is well visible in X-rays. Simple one-step system eliminates the need of using needle, guidewire, or drain placement radiological control. Safe, atraumatic drain prevents wound perforation. It is used for fluid drainage from the pleura, peritoneum, and cysts. Drains 9F and bigger allow for abscess drainage.

Reference:	Catheter size [F]	Catheter length [cm]	Catherer type	Two-pieces needle size [G]
OSD6F	6	26	PIGTAIL	18
OSD6FP			STRAIGHT	
OSD9F	9		PIGTAIL	15
OSD9FP			STRAIGHT	
OSD12FP	12		PIGTAIL	
OSD12F			STRAIGHT	
OSD14F	14		STRAIGHT	
OSD14FP			PIGTAIL	
OSD16F	16		STRAIGHT	
OSD16FP			PIGTAIL	

When ordering use Ref. Number, e.g.: OSD16F.

