

SPINAL ANAESTHESIA NEEDLE

Standard (QUINCKE point)



Set components:

- guiding needle
- spinal needle standard (Quincke point)

Reference:	Guiding needle size [G]	Spinal needle STANDARD size [G]	Spinal needle STANDARD length [mm]	Hub colour
IPPS18G50	---	18	50	pink
IPPS18G90	---		90	
IPPS19G50	---	19	50	cream
IPPS19G75	---		75	
IPPS19G90	---		90	
IPPS20G50	---	20	50	yellow
IPPS20G75	---		75	
IPPS20G90	---		90	
IPPS21G50	---	21	50	green
IPPS21G90	---		90	
IPPS22G50	---	22	50	black
IPPS22G75	---		75	
IPPS22G90	---		90	
IPPS23G90	---	23	90	blue
IPPS24G90	---	24	90	violet
IPPS25G50	20	25	50	orange
IPPS25G75			75	
IPPS25G90			90	
IPPS26G50	22	26	50	brown
IPPS26G90			90	
IPPS26G130			130	
IPPS27G90	22	27	90	grey
IPPS27G130		27	130	
IPPS29G90		29	90	

When ordering use Ref. Number, e.g.: IPPS27G90.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.

