

SPINAL ANAESTHESIA BLOCK SET

Standard type



Set components:

- STANDARD needle (Quincke type)
- guiding needle
- subcutaneous anaesthesia needle 0.5 x 25 mm
- needle for medicine drawing 0.9 x 40 mm
- regular syringe 2.5 ml
- regular syringe 5 ml

Reference:	Guiding needle size [G]	Spinal needle STANDARD size [G]	Spinal needle STANDARD length [mm]	Hub colour
ZIPPS18G	---	18	90	pink
ZIPPS19G	---	19		cream
ZIPPS20G	---	20		yellow
ZIPPS21G	---	21		green
ZIPPS22G	---	22		black
ZIPPS23G	---	23		blue
ZIPPS24G	20	24		blue
ZIPPS25G	20	25	90, 130	orange
ZIPPS26G		26		brown
ZIPPS27G	22	27		grey
ZIPPS29G		29	90	red

When ordering use Ref. Number, e.g.: ZIPPS27G.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.

