

EPIDURAL ANAESTHESIA SET SMALL

Set components:

- Tuohy needle
- epidural catheter with hub
- flat epidural filter 0.2 μm with self-adhesive element
- sleeve



Optional component:



system
for catheter fixation



10 ml low resistance syringe

The epidural set allows you to administer anesthetics continuously or repeatedly.

Reference:	Catheter size [G]	Tuohy needle size [G]	Tuohy needle length [mm]
ZZOM19G	22 (pediatric)	19	50
ZZOM18G	20	18	80
ZZOM17G	19	17	80
ZZOM16G	18	16	80

When ordering use Ref. Number, e.g.: ZZOM18G.

NOTE: The set can be purchased together with a low-resistance syringe.

While ordering, add the letters "SN" to the reference number – for example: ZZOM18GSN.

Special order can contain epidural catheter with SOFT tip (comes only with ZZOM18G).

While ordering, add the letter "S" to the reference number – for example: ZZOM18GS.

Special order can contain system for catheter fixation.

While ordering, add the letters "EM" to the reference number – for example: ZZOM18GEM.

