



## EMBOLECTOMY AND THROMBECTOMY CATHETERS

double lumen



The embolectomy catheters are used for the safe and effective removal of thromboembolic material from the peripheral arterial system. Double lumen catheter can be used to make the emboli visible with radiography by first inflating the balloon and then injecting radiopaque product through the distal orifice.

Reference:	Catheter size [F]	Catheter length [cm]	Diameter of filled balloon [mm]	Balloon volume filled with liquid [ml]
EFBD3F40	3	40	6	0.40
EFBD3F80		80		
EFBD4F40	4	40	8	0.50
EFBD4F80		80		
EFBD5F40	5	40	10	0.60
EFBD5F80		80		
EFBD6F40	6	40	12	0.80
EFBD6F80		80		
EFBD7F40	7	40	14	1.50
EFBD7F80		80		
EFBD8F40	8	40	16	3.00
EFBD8F80		80		

When ordering use Ref. Number, e.g.: EFBD7F80.

