

SPINAL ANAESTHESIA BLOCK SET

with Pencil-Point needle



Set components:

- PENCIL-POINT needle (Whitacre type)
- guiding needle
- subcutaneous anaesthesia needle 0.5 x 25 mm
- needle for medicine drawing 0.9 x 40 mm
- regular syringe 2.5 ml
- regular syringe 5 ml

Reference:	Guiding needle size [G]	Spinal needle PENCIL-POINT size [G]	Spinal needle PENCIL-POINT length [mm]	Hub colour
ZIPPW22G	- - -	22	90	black
ZIPPW24G	20	24	90	violet
ZIPPW25G	20	25	90	orange
ZIPPW25G120	20	25	120	orange
ZIPPW26G	20	26	90	brown
ZIPPW26G120	20	26	120	brown
ZIPPW27G	22	27	90	grey
ZIPPW27G120	22	27	120	grey

When ordering use Ref. Number, e.g.: ZIPPW27G.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.

