

EPIDURAL ANAESTHESIA SET ADVANCED



Set components:

- Tuohy needle
- epidural catheter with hub
- flat epidural filter 0.2 μm with self-adhesive element
- needle for drug administration 0.9 x 40 mm
- needle for subcutaneous anaesthesia 0.5 x 25 mm
- 10 ml low resistance syringe
- 10 ml regular syringe
- sleeve



Optional component:

system
for catheter fixation

The epidural set allows you to administer anesthetics continuously or repeatedly.

Reference:	Catheter size [G]	Tuohy needle size [G]	Tuohy needle length [mm]
ZZOR19G	22 (pediatric)	19	50
ZZOR18G	20	18	80
ZZOR17G	19	17	80
ZZOR16G	18	16	80

When ordering use Ref. Number, e.g.: ZZOR18G.

NOTE: Special order can contain epidural catheter with SOFT tip (comes only with ZZOR18G).

While ordering, add the letter "S" to the reference number – for example: ZZOR18GS.

Special order can contain system for catheter fixation.

While ordering, add the letters "EM" to the reference number – for example: ZZOR18GEM.

