

## DIALYSIS SET


single lumen, straight



### Set components:

- single lumen straight catheter
- needle 18G, length 7 cm
- guidewire
- dilator I
- dilator II

### Optional components



Reference: guidewire "J" type	Reference: guidewire Straight type	Catheter size [F]	Catheter length [cm]	Dilator I/II [F]/[F]	Guidewire size [inch]	Guidewire length [cm]
ZDJP8F10	ZDJP8F10P	8	10	8/10	0.035"	60
ZDJP8F15	ZDJP8F15P		15			
ZDJP8F20	ZDJP8F20P		20			

When ordering use Ref. Number, e.g.: ZDJP8F20.

#### NOTE:

On demand, it is possible to select additional components.

While ordering, add the following endings to the reference index:

**S** – scalpel

**Y** – Y hub

**SR** – Guiding syringe 5 ml

**SL** – Syringe 5 ml

**KLLT** – TEGO connector

**KLLI** – Injection cap

**IY** – V needle

**N** – Nitinol guidewire

**H** – Hydrophilic coating

E.g.: ZDJP8F15HYKLLT.

Different dialysis catheter sizes and lengths available on demand. For detailed information, please contact the company's representative.

