



EMBOLECTOMY AND THROMBECTOMY CATHETERS

guidewire-compatible double lumen



The embolectomy catheters are used for the safe and effective removal of thromboembolic material from the peripheral arterial system.

Reference:	Catheter size [F]	Catheter length [cm]	Diameter of filled balloon [mm]	Balloon volume filled with liquid [ml]	Use with guidewire [inch]
EFBDP3F40	3	40	6	0.40	0.016"
EFBDP3F80		80			
EFBDP4F40	4	40	8	0.50	0.018"
EFBDP4F80		80			
EFBDP5F40	5	40	10	0.60	0.025"
EFBDP5F80		80			
EFBDP6F40	6	40	12	0.80	0.027"
EFBDP6F80		80			
EFBDP7F40	7	40	14	1.50	0.035"
EFBDP7F80		80			
EFBDP8F40	8	40	16	3.00	0.038"
EFBDP8F80		80			

When ordering use Ref. Number, e.g.: EFBDP8F40.

