


## EMBOLECTOMY AND THROMBECTOMY CATHETERS

single lumen



The embolectomy catheters are used for the safe and effective removal of thromboembolic material from the peripheral arterial system.

| Reference: | Catheter size [F] | Catheter length [cm] | Hub colour | Diameter of filled balloon [mm] | Balloon volume filled with liquid [ml] |
|------------|-------------------|----------------------|------------|---------------------------------|--|
| EFB2F40    | 2                 | 40                   | white      | 5                               | 0.15                                   |
| EFB2F60    |                   | 60                   |            |                                 |  |
| EFB2F80    |                   | 80                   |            |                                 |  |
| EFB3F40    | 3                 | 40                   | green      | 6                               | 0.20                                   |
| EFB3F60    |                   | 60                   |            |                                 |  |
| EFB3F80    |                   | 80                   |            |                                 |  |
| EFB4F40    | 4                 | 40                   | red        | 8                               | 0.70                                   |
| EFB4F60    |                   | 60                   |            |                                 |  |
| EFB4F80    |                   | 80                   |            |                                 |  |
| EFB5F40    | 5                 | 40                   | white      | 10                              | 1.50                                   |
| EFB5F60    |                   | 60                   |            |                                 |  |
| EFB5F80    |                   | 80                   |            |                                 |  |
| EFB6F40    | 6                 | 40                   | blue       | 12                              | 2.00                                   |
| EFB6F80    |                   | 80                   |            |                                 |  |
| EFB7F40    | 7                 | 40                   | yellow     | 14                              | 2.50                                   |
| EFB7F80    |                   | 80                   |            |                                 |  |
| EFB8F40    | 8                 | 40                   | brown      | 16                              | 3.50                                   |
| EFB8F80    |                   | 80                   |            |                                 |  |
| EFB10F40   | 10                | 40                   | grey       | 20                              | 4.50                                   |
| EFB10F80   |                   | 80                   |            |                                 |  |

When ordering use Ref. Number, e.g.: EFB5F40.

Availability and specifications are subject to change. For more information, please contact your Balton sales representative.

